




2008 USBC TEAM USA TRIALS MEDIUM



DISTANCE:	39 FEET	VOLUME:	23.48 mL
RATIO:	2.14:1	FORWARD:	9.36 mL
DROP BRUSH:	39 FEET	REVERSE:	14.12 mL
TANKS:	KEGEL	PUMP:	40μL





FORWARD LOADS DATA

#	START	STOP	LOADS	MICS	SPEED	BUFF	TANK	DISTANCE	T.OIL
1	2L	2R	3	40	14	3	A - Kegel	0 → 4	4,440
2	6L	5R	1	40	14	3	A - Kegel	4 → 6	1,200
3	9L	7R	1	40	14	3	A - Kegel	6 → 8	1,000
4	11L	9R	1	40	18	3	A - Kegel	8 → 10	840
5	12L	10R	1	40	18	3	A - Kegel	10 → 13	760
6	14L	11R	1	40	18	3	A - Kegel	13 → 15	640
7	16L	13R	1	40	18	3	A - Kegel	15 → 18	480
8	2L	2R	0	40	18	3	A - Kegel	18 → 24	0
9	2L	2R	0	40	22	3	A - Kegel	24 → 31	0
10	2L	2R	0	40	30	3	A - Kegel	31 → 39	0

REVERSE LOADS DATA

#	START	STOP	LOADS	MICS	SPEED	BUFF	TANK	DISTANCE	T.OIL
1	2L	2R	0	40	30	3	A - Kegel	39 → 32	0
2	14L	12R	1	40	18	3	A - Kegel	32 → 29	600
3	12L	10R	1	40	18	3	A - Kegel	29 → 27	760
4	11L	9R	1	40	14	3	A - Kegel	27 → 25	840
5	9L	8R	1	40	14	3	A - Kegel	25 → 23	960
6	6L	6R	2	40	14	3	A - Kegel	23 → 19	2,320
7	2L	5R	2	40	10	3	A - Kegel	19 → 16	2,720
8	2L	2R	4	40	10	3	A - Kegel	16 → 11	5,920
9	2L	2R	0	40	10	3	A - Kegel	11 → 0	0



KEGEL PATTERN LIBRARY APP



PATTERNLIBRARY.KEGEL.NET



SCAN QR CODE TO
SHARE AND
DOWNLOAD
THIS PATTERN

